



Reference: PTXBC001201

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SIAT7-E; SIAT7E; ST6GalNAcV; alpha-N-acetylgalactosaminide alpha-2,6-sialyltransferase 5; GD1 alpha synthase; GalNAc alpha-2,6-sialyltransferase V; ST6 (alpha-N-acetyl-neuraminy-2,3-beta-galactosyl-1,3)-N-acetylgalactosaminide alpha-2,6-sialyltransferase		

For research use only.

GENE DESCRIPTION	ST6 (alpha-N-acetyl-neuraminyl-2,3-beta-galactosyl-1,3)-N-acetylglactosaminide alpha-2,6-sialyltransferase 5
GENE ID	81849
ORIGIN SPECIES	Human
ACCESSIONS	BC001201
VECTOR	pDONR223
NCBI SYMBOL	ST6GALNAC5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]