

R3HDM1-R3H domain containing 1 Gene

Reference: PTXBC001217

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	R3HDM; R3H domain-containing protein 1; R3H domain (binds single-stranded nucleic acids) containing; R3H domain containing 1		

For research use only.

GENE DESCRIPTION	R3H domain containing 1
GENE ID	23518
ORIGIN SPECIES	Human
ACCESSIONS	BC001217
VECTOR	pDONR223
NCBI SYMBOL	R3HDM1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggatgtctgatactgttactgtaaaagatgaaactgcaacaatgaaggattggaggcagaagtgaagataaccagagttgaaatcttatcaaatcagaaaactatgggaagatttggtagagaagaatgaacattgtattgagaacaatatagattgcaggtaaacaggaattttctctggctacaaataataattcagtatatatattgtttattttacatacaaaaactcacttctcaaaattctaaaagggttaagtga