



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

ZNF581-zinc finger protein 581 Gene

Reference: PTXBC001219

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

HSPC189; zinc finger protein 581

For research use only.

GENE DESCRIPTION	zinc finger protein 581
GENE ID	51545
ORIGIN SPECIES	Human
ACCESSIONS	BC001219
VECTOR	pDONR223
NCBI SYMBOL	ZNF581
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctggtgctgccatccccctgccctacagccttggcattttctcgtgagacataggaggccctccccctggacttgcgctccccgaacctggacctctctctccatcggaatctccccaggctcatctctcccaaggcccaacctactcgttattgaaactcagggtgtccc
ctacacagtgctggfggagcaggatcacagaggaggcggggccagtggggctcaggccgagaaaagtgtacagctgcgccctgtgtctcaagggtcttcgagtacattctcattctcagcgacacagatcacaccactcggaggtaagcccttcgagtgtgacatctgtggg
aaggcttataagcgcgcagccacttggcacggccacattcattcacctgggtgtgtggcgcccccagctgccgcctctgccttcgcgccttgggggtgagctgtgccacgacagccgggtgcactctggggagcgccctgactgttccacactgcctctgc
cgctttatggagcagaacacatcagaaaacacgcggatgtggaagcatcatgga