

APH1A-anterior pharynx defective 1 homolog A (C. elegans) Gene

Reference: PTXBC001230

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	APH1A gamma secretase subunit; 6530402N02Rik; APH-1; APH-1A; CGI-78; gamma-secretase subunit APH-1A; anterior pharynx defective 1 homolog A; aph-1alpha; presenilin-stabilization factor; aph-1 homolog A, gamma-secretase subunit		

For research use only.

GENE DESCRIPTION	anterior pharynx defective 1 homolog A (C. elegans)
GENE ID	51107
ORIGIN SPECIES	Human
ACCESSIONS	BC001230
VECTOR	pDONR223
NCBI SYMBOL	APH1A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggggctgcgggtttttcggctgcacttctcgtcgttcgcccggccttcgcgtttcttgatcactgtggctggggaccgcttcggtatcatcctggtcgcaggggcattttctggctggtctccctgctcctggcctctgtggtctggttcattctgtccatgtgaccgaccgtca
gatgccgggctccagtacggcctcctgatttttggtgctgctgtctctctacaggaggtgttcgcttgcctactacaagctgcttaagaaggcagatgaggggttagcatcgtgagtgaggacggaagatcacccatctccatccgagatggcctatgtttctggtctctccttcg
gtatcatcagtggtgtctctctgttatcaatattttggctgatgcacttggccagggtgtggttgggatccatggagactcacccattacttctgacttcagcctttctgacagcagcattatcctgctccatccttttggggagttgttctttgatgcctgtgagaggagacgggtactggcctt
tgggcctggtgttgggagtcacactactgacatcgggactgacattctgaacccctggtatgaggccagcctgctgccatctatgcagtcactgtttccatggggctctgggccttcacacagctggagggtccatccgaagtattcagcgcagcctctgtgccgacggcaggaggac
agtcgggtgatggttatcttgcctgcgcacccaccgaggactga