

GSTK1-glutathione S-transferase kappa 1 Gene

Reference: PTXBC001231

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GSTK1-1; GST; GST 13-13; GST13; GST13-13; hGSTK1; glutathione S-transferase kappa 1; GST class-kappa; glutathione S-transferase k1; glutathione S-transferase subunit 13 homolog		

For research use only.

GENE DESCRIPTION	glutathione S-transferase kappa 1
GENE ID	373156
ORIGIN SPECIES	Human
ACCESSIONS	BC001231
VECTOR	pDONR223
NCBI SYMBOL	GSTK1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggccctgcccgcaccgtggagctcttctatgacgtgctgtccccctactcctggctggcttcgagatcctgtgccggtatcagaatatctggaacatcaacctgcagttgcggccagcctcataacaggatcatgaaagacagtggaaacaagcctcaggtctgttccccg
caaggactatacatggcaatgacttaagctcctgagacaccatccagattccatccacttcccgaagatttctgtctgtgctgtgcttgaagaaaggagttgtctgcatgctgttctcaccgctgaaactggagcatccagagatgctggagaaagcgtccgggagctgtgga
tgcgctgtgtgcaaggaatgaagacatcaccgagcgcagagcatcctggcggctgcagagaaggctggtatgtctgcagaacaagcccaggacttctgaaaagatcgaacgccaagggtgaagaaccagctcaaggagaccactgaggcagcctgcagataggagcctt
gggctgccatcaccgtggccatgtggatggccaaaccacatgttatttggcttgaccggatggagctgctggcgacacgtctgggagagaagtggatggccctatactccagccgtgaatgccagactttaa