

ARHGAP17-Rho GTPase activating protein 17 Gene

Reference: PTXBC001241

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MST066; MST110; MSTP038; MSTP066; MSTP110; NADRIN; PP367; PP4534; RICH-1; RICH1; WBP15; rho GTPase-activating protein 17; RhoGAP interacting with CIP4 homologs 1; neuron-associated developmentally regulated protein; rho-type GTPase-activating protein 17;		

For research use only.

GENE DESCRIPTION	Rho GTPase activating protein 17
GENE ID	55114
ORIGIN SPECIES	Human
ACCESSIONS	BC001241
VECTOR	pDONR223
NCBI SYMBOL	ARHGAP17
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcacccaacccaatagccaggccctcccaaccccatggcattgccagtgagcatggacttgagcagccatctcacacccctcccagactccaagcgcgccctagctactcgcccttaggaaaacagaaccccgactgtgccagctcctcagaccctggcaggggtaaccctg
aaactgcacagccacatgctggaaccttaccgagaccgagaccagtaccaaagccaaggaaccggccagcgtgccccaccccccaacctctgtgtccactcagctggggacagcagctcaccaacacagcaccacacagctccaagatagtaacagatgtatga