



Reference: PTXBC001261

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM50B; D6S2654E; XAP5-like protein; family with sequence similarity 50 member B		

For research use only.

GENE DESCRIPTION	family with sequence similarity 50, member B
GENE ID	26240
ORIGIN SPECIES	Human
ACCESSIONS	BC001261
VECTOR	pDONR223
NCBI SYMBOL	FAM50B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1