

NMI-N-myc (and STAT) interactor Gene

Reference: PTXBC001268

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	N-myc-interactor; N-myc and STAT interactor		

For research use only.

GENE DESCRIPTION	N-myc (and STAT) interactor
GENE ID	9111
ORIGIN SPECIES	Human
ACCESSIONS	BC001268
VECTOR	pDONR223
NCBI SYMBOL	NMI
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaagctgataaagatgacacacaacaattcttaaggagcattcgccagatgaatttataaaagatgaacaaataagggaactaattgatgaattacaagaataatcaactaaagaaggagatccaaagcttgaacggaggtacaagaggtaccagaattccagattaaa
gaggatattcctgaacaaagatgaattcttcatcgttgaaactcctgagaatgacagccagttgtcaaatatcctctgttcattcaagtgagctcgaaagttccttatgagatacaaaaaggacacacttatcactttgaaaaagaagaagttgctcaaatgtggaagcatgagtaaac
atcatgtacagataaaagatgtaaatctggaggttacggccaagccagttccataaattcaggagtcagattccaggtttatgtagaagtttctaaaatgaaaatcaatgttactgaaattcctgacacattgcgtgaagatcaaatgagagacaaactagagctgagcttttcaagtcctgaaa
tggaggcggagaggtggaccgctggactatgacagacagtcgggagtcagtcacagtttgggagattggagtggtgacaagattttgaaaaagaagaataccctctttataataatcaaacctgccaatagattctgtttctccatcacagaatacacttgaaaaagtatcag
atatcttcaggaaacatctaagaggacagtgctctgacaggaatggaaggcattcaaatggatgaagaattgtggaggatttaataacattcactttcaacgggcaagaatggaggtggagaagtagatgtgtcaagtgtctctaggtcaacctacatagcatactttgaagaatag