

## CDC2L5-cell division cycle 2-like 5 (cholinesterase-related cell division controller) Gene

Reference: PTXBC001274

|                |                                                                                                                                                                                                                                                   |                     |                |
|----------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                                               | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                                              | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | CDC2L5; CDC2L; CHED; hCDK13; cyclin-dependent kinase 13; CDC2-related protein kinase 5; cell division cycle 2-like protein kinase 5; cell division protein kinase 13; cholinesterase-related cell division controller; cyclin dependent kinase 13 |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | cell division cycle 2-like 5 (cholinesterase-related cell division controller)                    |
| <b>GENE ID</b>          | 8621                                                                                              |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC001274                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | CDC2L5                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgctgcctgaagataaagaagctgatacgtacgaggaatatttcagtaaaagcagttaaaaaagaagtagaaaaaagcagctgctcttctgctgattaccgctgccccctgagctaccaggaggagatgatcttcaagagtcagaggaaaaaagacacacagttacat  
agtaaaaggaggcctaaaatatgtggcctcgctatggtgaaccaaagaaaaagatattgactggggaaaaacgctgctggaataattgataatcaggaattattggagaagggtactacggacaagtttacaagccagggaataaagacactggagaaatggtagccttaaaaaaagt  
acgtctggataatgaaaggaaggccttccaattacagcaattcagagaaattaaaattctcggcagcttaccatcagagattatcaatatgaaggaatagtgactgataaagaagatgctttggattcaagaaggacaaaggtgcatfctgtgttgaatataatggacatgatctga  
tgggactactggaatcaggcttgggttcattttaaataaaaatcacataaagtcatttatgagacagctcatggaggctggtatttgcataagaagaacttttgcataagataatgaatgttccaatccttctaataatagagggcagataaaaacttcagactttggacttgcctgattgtat  
agctcagaagaagtcggcgtataactaacaaggtaattactttatggtaccgtccacgtgaactgctactgggagaagaacgatacacaccagcattgatgatggagctgtggctgtatccttggcgaactcttcaataaaacctatatttcaagcaaatcagggaacttgcacaactaga  
attaataagacatgaagaaaatgaagttctgataaacaactga