

HMGN4-high mobility group nucleosomal binding domain 4 Gene

Reference: PTXBC001282

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HMG17L3; NHC; high mobility group nucleosome-binding domain-containing protein 4; high mobility group protein N4; high-mobility group (nonhistone chromosomal) protein 17-like 3; high-mobility group protein 17-like 3; non-histone chromosomal protein HMG-17		

For research use only.

GENE DESCRIPTION	high mobility group nucleosomal binding domain 4
GENE ID	10473
ORIGIN SPECIES	Human
ACCESSIONS	BC001282
VECTOR	pDONR223
NCBI SYMBOL	HMGN4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcccaagagaaaggcaaaaggagatgctaaagggtataaagcaaaagggtgaaggatgagccacagaggagatcagctcggttgctgctaaaccagctctccaaaccagagccaggcctaaaaaggcctctgcaaaaggaggagagaagcttccaaagggagaaaggga
aagcagatgctggaaaggatggaaacaacctgcaaaaaaccgagatgcctctacactccagtcctccagaaagcggaaggcactggggatgccaaagtga