

C2orf49-chromosome 2 open reading frame 49 Gene

Reference: PTXBC001310

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	asw; ashwin; chromosome 2 open reading frame 49		
For research use only.			

GENE DESCRIPTION	chromosome 2 open reading frame 49
GENE ID	79074
ORIGIN SPECIES	Human
ACCESSIONS	BC001310
VECTOR	pDONR223
NCBI SYMBOL	C2orf49
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgggggatgtggcggtcgagctgcacggactcgaactgctgctgacccggagctgctgccaggagttcctctctcaactctggagcagaagaacatagctgttgaactgatgaagagtaacaagacagctctactgacctttatgtccaacatgcaataccattgcct
cagagggaatttccgaagaatagatggggaaaatgatgaaaagaaaagagaacaacatgagattaaaaatgagactaaaaaggagtagcactgtatgggttaaggaaaagaccctcatcgtattgatggaagtcaacaagtacaagcataaaagtgaagaaacagagagaatgg
agataatgatcgactgaagctccccgcagcgaagctttaccagtaattgctttgaaaaattatcaaatctcttcgagtggtttacccttaatttgcctccaatttgcctgtgaacaataaaacggaacacaataataatgacgctaacaagacacatgactaacgcataaggaagtcctt
caggccctgtgaagtcgccaccattgtccctgttggaactactccagtgaagttaagagagctgctcctaagaagaggcagagccatgaataacctgaagccccacaagcaaaagggaagataacaacatgttacttggccctga