

39699-septin 8 Gene

Reference: PTXBC001329

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	septin 8; SEP2; septin-8		

For research use only.

GENE DESCRIPTION	septin 8
GENE ID	23176
ORIGIN SPECIES	Human
ACCESSIONS	BC001329
VECTOR	pDONR223
NCBI SYMBOL	39699
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagaactagacagcaaggtgaacattatcccatcatcgccaaggctgacacatctccaagagcgagctccacaagtcaagatcaagatcatggcgaggttggtcagcaacgggtccagatctaccagttcccccggatgatgaggtgttgcagagattaacgcagtcatga
atgcacatctgcccttggcgttggtggcagcaccgaggaggtgaaggtgggaacaagctgggtccgagcacggcagtlacccctggggagtggtgaggtggagaatgagaatcactgcgactcgtgaagctgcgggagatgttgatccgggtgaacatggaagacctccgcgag
cagacccacagccggcactacgagctctaccgctgcaagttggaggagatgggtttcaggacagcgtggtgacagccagcccttcagcctacaagagacatacaggccaagaggaaggagttcctaagtgcgctgcagaggaaggaggaagagatgaggcagatgtttgt
caacaaagtgaaggagacagagctggagctgaaggagaaggaaaggagctccatgagaagtttgagcacctgaagcgggtccaccaggaggagaagcgcaaggtgaggaaaagcgccgggaactggaggaggagaccaacgccttcaatgccggaaggctgcgtgga
ggcctgcagtcgagggccttgcagccacctgcagcagccctgagggaaggacaaggacaagaagaattaa