

POLD4-polymerase (DNA-directed), delta 4 Gene

Reference: PTXBC001334

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	POLDS; p12; DNA polymerase delta subunit 4; DNA polymerase delta smallest subunit p12; polymerase (DNA) delta 4, accessory subunit; polymerase (DNA-directed), delta 4, accessory subunit; DNA polymerase delta 4, accessory subunit		

For research use only.

GENE DESCRIPTION	polymerase (DNA-directed), delta 4
GENE ID	57804
ORIGIN SPECIES	Human
ACCESSIONS	BC001334
VECTOR	pDONR223
NCBI SYMBOL	POLD4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccggaagcggtctcactgattcctaccggttggaagagaggaggccgctgggcacagcaaggggagctggcaccgagctaggtctctga