

LRRC61-leucine rich repeat containing 61 Gene

Reference: PTXBC001354

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSPC295; leucine-rich repeat-containing protein 61; leucine rich repeat containing 61		
For research use only.			

GENE DESCRIPTION	leucine rich repeat containing 61
GENE ID	65999
ORIGIN SPECIES	Human
ACCESSIONS	BC001354
VECTOR	pDONR223
NCBI SYMBOL	LRRC61
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaccctccagcggagaagccggagaggtggcggactgcagatcacacccagctgctgaagtcacgcacggcaggtctccctggagtcacatcctgctactgaagctgcgtggcttgggactggctgacctgggctgcctgggagagtgctgggctggagtgctggac
ctatcaggcaacgcgtcacccacctgggcccgtggcctcctfgcgcagctagctgtgtcaatgtctccaacaatcggtgacggcctggagccactggccacctgfgagaacttgcagagttctcaatgccgcaggcaacctactggccacccgggcccagctgcagtgctggc
tgggctaccgtgctggagtaactgcgtccgagaccctfggccggctcagcaaccgctctgtgccaacccctcctactgggctgcagtcgggagctgctgctggcctgaaagtcacgacgtgagcgtgtgattggcgtgtgtagtgaattctaccagctgtgccgagacct
ggacagctcctgtgctccagctccagtcaggccccagagccaccgaggccagccctgggtggagccaggctactgggagtcctggccagccggagcagctccatcctggaggaggcctgccggcagttccaggacacactgcaggagtgctgggacctggaccgagggc
cagcgacagctggccagggcggagcaggtactcagctctgcggccacacctctctctctga