

ZNF747-zinc finger protein 747 Gene

Reference: PTXBC001361

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KRAB domain-containing protein ZNF747; zinc finger protein 747		
For research use only.			

GENE DESCRIPTION	zinc finger protein 747
GENE ID	65988
ORIGIN SPECIES	Human
ACCESSIONS	BC001361
VECTOR	pDONR223
NCBI SYMBOL	ZNF747
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaccgatccgtcgtggggctgacagtcacatggcgccgtctggcccgctccctcccggaacccaacggggggatccgagtgagaaaagccgggcccgtgagcttcgcccagctggccgtgtactctccgggaggagtgaggctgctgcggccgcgcaga
gggcccgtgacgggacgtgatggggagacacacggccacctggcgccgtcgtgagagcccacctgctgctgggcccctgcgcctccacagccctgcccgcctctgggagctgcgtgtgagtggtggggcccccggggccggcagggcgcctcctgcagcgtgggg
tttcggttcttccccagggagtcggaggcagcaagccggcgctcatctcctgggtggaggagaaagccgaactgtgggatccggctgcccaggatccggaggtggcgaagtgtccgacagaagcggaccagcagattccagaacaaggaggaagacaaggggaagg
gacgggagccctggagaagcccacccctgtggcccggtctcctgggctga