

SEC22B-SEC22 vesicle trafficking protein homolog B (*S. cerevisiae*) Gene

Reference: PTXBC001364

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	vesicle-trafficking protein SEC22b; ERS-24; SEC22L1; ER-Golgi SNARE of 24 kDa; SEC22 vesicle trafficking protein homolog B; SEC22 vesicle trafficking protein-like 1; SEC22 homolog B, vesicle trafficking protein (gene/pseudogene)		

For research use only.

GENE DESCRIPTION	SEC22 vesicle trafficking protein homolog B (<i>S. cerevisiae</i>)
GENE ID	9554
ORIGIN SPECIES	Human
ACCESSIONS	BC001364
VECTOR	pDONR223
NCBI SYMBOL	SEC22B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgttgctacaatgatcgcccgagtgccggacgggctcccgctggccgctcgcagcaggaggacgaacagctctggccggacctcaacagatcagagtcaggcctaagcaactcttcgaaagttgaatgaacagtcctccaccagatgtaccttgaagcaggagccatgact
tttccactacattatgagcaggggggtgtgtatttggttttatgtgaagctgcctccctaagaagttggcttttgcctaccatagaagatttgcactcagaatttgatgaacagcatggaaagaagggtgccactgtgtcccgacctatctcttattgaatttgatactttcaatcagaaaaccaagaag
ctctacattgacagctgtgtcgaagaatctaggctccatcaacactgaattgcaagatgtgcagaggatcatgtgtggccaatttgaagaaggtttacaacgaggagaagcactctcagcattggattcaaaaggctaaacttgcagctgtgtccaagaataccgccaggtgcgaag
tactgaacatgcgttccactatgccaactgcagcagtagctgtatttttcatcatgttaatagtgtatgtccgattctgtgctgtga