

C9orf95-chromosome 9 open reading frame 95 Gene

Reference: PTXBC001366

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C9orf95; NRK1; bA235O14.2; nicotinamide riboside kinase 1; NRK 1; RNK 1; nicotinic acid riboside kinase 1; nmR-K 1; ribosylnicotinamide kinase 1; ribosylnicotinic acid kinase 1		

For research use only.

GENE DESCRIPTION	chromosome 9 open reading frame 95
GENE ID	54981
ORIGIN SPECIES	Human
ACCESSIONS	BC001366
VECTOR	pDONR223
NCBI SYMBOL	C9orf95
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaaacatttatcattggaatcagtggtgtgacaaacagtggcaaaacaacactggcctaagaattgcagaacacaccccccaattgcagtgatcatatcaggatgattcttcaagccagagctgagatagagacagataaaatggattttgcagtacgatgtgcttgaagcacttaac
atggaaaaaafgatgcagccatttctgctggatggaaagcgcaagacactctgtgtatcaacagaccaggaagtgctgaggaaattccatttaacatcgaaggtttcttcttttaattataagcccttgacactatatggaatagaagctatttctgacgattccatatgaagaatg
aaaaggaggaggagtacaagggtctatcagccctccagactctccgggatactttgatggccatgtgtggccatgtatctaaagtacagacaagaatgcaggacatcacatgggaagttgtgtacctggatggaacaaaatctgaagaggacctcttttgcagtatatgaagatctaata
caagaactagcaaaagtggttgcagtgacagcataa