

SNRNP25-small nuclear ribonucleoprotein 25kDa (U11/U12) Gene

Reference: PTXBC001381

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C16orf33; U11/U12 small nuclear ribonucleoprotein 25 kDa protein; U11/U12 snRNP 25 kDa protein; U11/U12 snRNP 25K; minus-99 protein; small nuclear ribonucleoprotein 25kDa (U11/U12); small nuclear ribonucleoprotein, U11/U12 25kDa subunit; small nuclear rib		

For research use only.

GENE DESCRIPTION	small nuclear ribonucleoprotein 25kDa (U11/U12)
GENE ID	79622
ORIGIN SPECIES	Human
ACCESSIONS	BC001381
VECTOR	pDONR223
NCBI SYMBOL	SNRNP25
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacgtgttcaggagggtctggtatggtggtcaggaccgctgctcgcgatctgccgatccaggttactctggaagaagtcactcccaaatagccctagatacggccagcgcaatgacggtccgagtgcaagatggatggagaagtaatgccctggttgtagtcagagt
gccacagtcctggacctgaagaaggccatccagagatagtgacgtcctcaagcaggagcgtgaagggggcattcagcacatcagctggtcctactgtgtggaggacgtaccatctgacctctgcaggagagaaactcacggaagacagaaagagctccgagactacggcatccgga
tcgagacgaggtttcctcatcaaaaagctgaggcaaaagtga