

UQCRQ-ubiquinol-cytochrome c reductase, complex III subunit VII, 9.5kDa Gene

Reference: PTXBC001390

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MC3DN4; QCR8; QP-C; QPC; UQCR7; cytochrome b-c1 complex subunit 8; complex III subunit 8; complex III subunit VIII; low molecular mass ubiquinone-binding protein (9.5kD); ubiquinol-cytochrome c reductase complex 9.5 kDa protein; ubiquinol-cytochrome c red		

For research use only.

GENE DESCRIPTION	ubiquinol-cytochrome c reductase, complex III subunit VII, 9.5kDa
GENE ID	27089
ORIGIN SPECIES	Human
ACCESSIONS	BC001390
VECTOR	pDONR223
NCBI SYMBOL	UQCRQ
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccgcgagtttgggaatctgacgcggatgcggcatgtatcagctacagctgtcaccgttcgagcagcgcgcctatccgacgtcttactaaaggaatcccaatgtctgcgccgattcgggagctcttcttcgcgttggtgccgagttttagtgttttatctatcacatggg
ggactgaagagttcgagagatccaagagggaagaatccagctgcctatgaaatgacaaatga