



Reference: PTXBC001395

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	WD40 protein Ciaol1; CIA1; WDR39; WD repeat domain 39; WD repeat-containing protein 39; cytosolic iron-sulfur protein assembly 1 homolog; cytosolic iron-sulfur assembly component 1		

For research use only.

GENE DESCRIPTION	cytosolic iron-sulfur protein assembly 1 homolog (S. cerevisiae)
GENE ID	9391
ORIGIN SPECIES	Human
ACCESSIONS	BC001395
VECTOR	pDONR223
NCBI SYMBOL	CIAO1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]