

S100A11-S100 calcium binding protein A11 Gene

Reference: PTXBC001410

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HEL-S-43; MLN70; S100C; protein S100-A11; MLN 70; calgizzarin; epididymis secretory protein Li 43; metastatic lymph node gene 70 protein; protein S100-C; S100 calcium binding protein A11		

For research use only.

GENE DESCRIPTION	S100 calcium binding protein A11
GENE ID	6282
ORIGIN SPECIES	Human
ACCESSIONS	BC001410
VECTOR	pDONR223
NCBI SYMBOL	S100A11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcaaaatctccagccctacagagactgagcgggtgcacgagtcctgattgctgtctccagaagtatgctggaaggatgggtataactacactctctccaagacagagttcctaagcttcatgaatacagaactagctgccttcacaagaaccagaagacccctgggtgccttgacc
 gcatgatgaagaactggacaccaacagtgatggcagctagattctcagaatttcttaactctgattggtggcctagctatggcttgccatgactccttctcaaggctgtcccttcccagaagcgggacctga