

TMEM54-transmembrane protein 54 Gene

Reference: PTXBC001418

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BCLP; CAC-1; CAC1; transmembrane protein 54; beta-casein-like protein; protein CAC-1		

For research use only.

GENE DESCRIPTION	transmembrane protein 54
GENE ID	113452
ORIGIN SPECIES	Human
ACCESSIONS	BC001418
VECTOR	pDONR223
NCBI SYMBOL	TMEM54
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgtctgcgcctggaggcctgagtggtggcgacttcggaaggtgctgatgaagacagccctggtgctggtgctggtggccatgtgagcttcacacagctgccctgttccatggcagagtgctgcgtctacgtgggcacccctcaagatgcgggtgctctgcagtactgctgtgtca
acatcctctctgtcactccgcatcgtggtcatcacttcaggcatcgagccatcgtgtgttcacgctacccctcctagcaccctcctgcgcctggacagtggttttagctcgagcgtggcctgtgctctcctttctctgacctgtgcccctggcctcttggcctccatgcctatgacctttgccacca
gggcaaggcactgctgctgctgctgctgttgggagctctgaactactggccctcgacactgactgctccttcgacccacacgacatttatagctccagcctgtgctctggggcctcgccttagtctctgctgtggcggaacgtgtttgctgtacgtgtgctcagctcaccaccagct
gctggagctgagccctggtgggggaaagcagccaccatgatgcgggagaaccagagctggtggaggccgtgacctgctgagctgcaccagctctgagcctctgacctctga