

UQCRH-ubiquinol-cytochrome c reductase hinge protein Gene

Reference: PTXBC001426

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	QCR6; UQCR8; cytochrome b-c1 complex subunit 6, mitochondrial; complex III subunit 6; cytochrome c1 non-heme 11 kDa protein; mitochondrial hinge protein; ubiquinol-cytochrome c reductase complex 11 kDa protein; ubiquinol-cytochrome c reductase, complex II		

For research use only.

GENE DESCRIPTION	ubiquinol-cytochrome c reductase hinge protein
GENE ID	7388
ORIGIN SPECIES	Human
ACCESSIONS	BC001426
VECTOR	pDONR223
NCBI SYMBOL	UQCRH
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggactggaggacgagcaaaagatgcttaccgaatccggagatcctgaggaggaggaagaggaggaggaattagtgatccctaacaacagtgagagagcaatcgagcagttggagaatgtgtaaaggccgggagcgctagagctctgtgatgagcgtgtatcctctcgatcacatacagaaggattgcacggaggagctcttgactcttgcatgacgagggaccattgcgtggcccaaaactctttaacaactgaaataa