



Reference: PTXBC001427

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CDG2C; FUCT1; GDP-fucose transporter 1; solute carrier family 35 (GDP-fucose transporter), member C1; solute carrier family 35 member C1		

For research use only.

GENE DESCRIPTION	solute carrier family 35, member C1
GENE ID	55343
ORIGIN SPECIES	Human
ACCESSIONS	BC001427
VECTOR	pDONR223
NCBI SYMBOL	SLC35C1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]