

POP7-processing of precursor 7, ribonuclease P/MRP subunit (*S. cerevisiae*) Gene

Reference: PTXBC001430

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	POP7 homolog, ribonuclease P/MRP subunit; ribonucleases P/MRP protein subunit POP7 homolog; POP7 (processing of precursor, <i>S. cerevisiae</i>) homolog; 0610037N12Rik; RPP2; RPP20; ribonuclease P protein subunit p20; RNaseP protein p20; hPOP7; processing of pre		

For research use only.

GENE DESCRIPTION	processing of precursor 7, ribonuclease P/MRP subunit (<i>S. cerevisiae</i>)
GENE ID	10248
ORIGIN SPECIES	Human
ACCESSIONS	BC001430
VECTOR	pDONR223
NCBI SYMBOL	POP7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagaaaaccgagagccccgcggtgctgtggaggctgaactggatccagtggaatacaccccttaggaaaaggcttccagccgctgccccggagaccgaatgacattatgtcaacatgaagacggacttaagccccagctggcccgctgccagaagctgctggacggaggg
gccccgggtcagaacgctgctgtgagatctacatcacggcttgggcttgccatcaaccgcccatacaatcgcgctgcagctgcaggcggcgagcttcgggtccttgaggtggtgccaatatccctccaccgtggagcttggatgagctggagccagagaccgacacacggg
agccactgactcgatccgcaacaactcagccatccatccgagcttcagggtcacaccaagtaa