

S100A6-S100 calcium binding protein A6 Gene

Reference: PTXBC001431

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	2A9; 5B10; CABP; CACY; PRA; protein S100-A6; MLN 4; calcyclin; growth factor-inducible protein 2A9; prolactin receptor-associated protein; S100 calcium binding protein A6		

For research use only.

GENE DESCRIPTION	S100 calcium binding protein A6
GENE ID	6277
ORIGIN SPECIES	Human
ACCESSIONS	BC001431
VECTOR	pDONR223
NCBI SYMBOL	S100A6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcatgccctggtcagggcattggcctcctcgtggccatcttcacaagtactccggcagggagggtgacaagcacaccctgagcaagaaggagctgaaggagctgatccagaaggagctcaccattggctcgaagctgcaggatgctgaaattgcaaggctgatggaagact
 tggaccggacaagaccaggaggtgaactccaggagtagtgcaccttctggggccttggccttgatctacaatgaagccctcaagggtga