

RBX1-ring-box 1 Gene

Reference: PTXBC001466

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	E3 ubiquitin-protein ligase RBX1; BA554C12.1; RNF75; ROC1; RING finger protein 75; RING-box protein 1; ZYP protein; regulator of cullins 1; ring-box 1, E3 ubiquitin protein ligase; ring-box 1		

For research use only.

GENE DESCRIPTION	ring-box 1
GENE ID	9978
ORIGIN SPECIES	Human
ACCESSIONS	BC001466
VECTOR	pDONR223
NCBI SYMBOL	RBX1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgagcagcgatggatgtggatacccgagcggcaccacacggcgcgccgcaagaagcgctttgaagtgaagaaagtggaatgcagtagccctctggcctgggatattgtggtgataactgtgccatctgcaggaaccacattatggatctttgcatagaatgtcaagctaaccag
gcgtccgctacttcagaagagtgtactgtcgcatgggaggtctgtaacctgcttttccactccactgcactctcgtgctcaaaacacgacaggtgtgtccattggacaacagagagtggaattccaaaagtagtggcactag