



Reference: PTXBC001471

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GNT1.2; GnT I.2; LGMD2O; MEB; MGAT1.2; RP76; gnT-I.2; protein O-linked-mannose beta-1,2-N-acetylglucosaminyltransferase 1; UDP-GlcNAc:alpha-D-mannoside beta-1,2-N-acetylglucosaminyltransferase I.2; protein O-linked mannose N-acetylglucosaminyltransferase		

For research use only.

GENE DESCRIPTION	protein O-linked mannose beta1,2-N-acetylglucosaminyltransferase
GENE ID	55624
ORIGIN SPECIES	Human
ACCESSIONS	BC001471
VECTOR	pDONR223
NCBI SYMBOL	POMGNT1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]