

TMED10-transmembrane emp24-like trafficking protein 10 (yeast) Gene

Reference: PTXBC001496

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	P24(DELTA); S31I125; S31III125; TMP21; Tmp-21-I; p23; p24d1; transmembrane emp24 domain-containing protein 10; 21 kDa transmembrane trafficking protein; p24 family protein delta-1; p24delta; p24delta1; testicular tissue protein Li 206; transmembrane emp24		

For research use only.

GENE DESCRIPTION	transmembrane emp24-like trafficking protein 10 (yeast)
GENE ID	10972
ORIGIN SPECIES	Human
ACCESSIONS	BC001496
VECTOR	pDONR223
NCBI SYMBOL	TMED10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctggttctgctgccaccagcccgccgcccctttccgttagcgttctgtctgttctgtctgtcccccagattggtcttgcctctcctccatctgccattactctgcaagtgcctccgtgaggagattcacaaggacctgctagtactggcgctacgagatctccgacc
agctctggggcgctggcgccctgcccagccacctcaagatcacagattctgtggccattctctactccaaaggagatgcaaccaaggggaaattgctttaccactgaagattatgacatgtttgaagtgtgtttgagagcaagggaacaggcgccgatactgaccaactctgatcc
tagacatgaagcatggagtggagcgaaaaattacgaagagattgcaaaagtgaagctcaaacattagaggtagagctgcgacgctagaagaccttcagaatctattgttaattgttcctacatgaagaagagagaagagagatgcgtgataccaacgagtcacaacacac
tcgggtctatactcagcatctttcaatgttctgtctcattggactagctacctggcaggtcttctacctgcgacgttctcaaggccaagaattgattgagtaa