

## TGFB1I1-transforming growth factor beta 1 induced transcript 1 Gene

Reference: PTXBC001507

|                |                                                                                                                                                                                                                                                                 |                     |                |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                                                             | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                                                            | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | ARA55; HIC5; TSC-5; transforming growth factor beta-1-induced transcript 1 protein; androgen receptor coactivator 55 kDa protein; androgen receptor coactivator ARA55; androgen receptor-associated protein of 55 kDa; hydrogen peroxide-inducible clone 5 prot |                     |                |

For research use only.

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | transforming growth factor beta 1 induced transcript 1                                           |
| <b>GENE ID</b>          | 7041                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC001507                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | TGFB1I1                                                                                          |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgccaaaggtcagggtcccaagagcgcctcgaggagcctctaccctccccatctatggccaccagccacagacaggggtctggggaggtctcaggagcctcgggggacaaggaccacctgtacagcacgggtatgaagcctcggtccccaagcctgcagccccggcgg  
ccccctccatctctctccagcgggtgtcttgggtaccgggtctgtgagctagatcgggtgtcttcagggaacttaatgccactcagttcaacatcacagatgaaatcgtctcagttcccatctagcaaggtggtctcaggagagcagaaggaggaccagctctgaagataagaaaagaccag  
cctccctccagcccgctctctggcctccaaagcctctgccacctcagccactctggagctggatagactgatggcctcactctctgactccgcgttcaaaacactctccagcctctgggccaactcagccaccgggtggtgagctccacaatgagggtccccatccccaccagac  
cgactggcaaggcagcctagacacatgctggggtgctgagctccgacctagccgccgggggtgttcccaccagggcgaaggcctctgtggtcctctgcaataaacctattgctgggcaagtgtgacggctctggcgccgcctggcaccggagcactctgttgcggaggtg  
ttccaccgcccctgggagggcagcagcttcttgagaaggtgagagccctctgccccgagtgctacttggagcgttctcgcgaagatgtggtctctgcaaccagccatccgacacaagatgtgacgccttgggcaactcactggcaccagagcatttctgctgctcagttgcggg  
gagcccttcggagatgagggtttccacgagcgcgaggccgccccctactgcccgggacttctgcagctgttgcgcccgctgccagggtgccaggcccatcctggataactacatctcggcgctcagcgcgctctggcaccgggactgttctgctgcagggaatgcttcg  
gcccttctcgggagggcagcttttgcagcacgaggccgcccgttctgagagaaccacttccagcagcagcggtctgctgtgcggcagctgtggcctcctgtgaccggccgctgcgtgtcggccctgggtcgccttccaccggaccacttccatgcacttctgctgcgc  
gctcaccaggggtctctccaggagcgcggccgcaagccctactgccagccctgcttctgaagctctcggctga