



Reference: PTXBC001513

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	VPS18, CORVET/HOPS core subunit; PEP3; vacuolar protein sorting-associated protein 18 homolog; hVPS18; vacuolar protein sorting 18 homolog; vacuolar protein sorting protein 18		

For research use only.

GENE DESCRIPTION	vacuolar protein sorting 18 homolog (S. cerevisiae)
GENE ID	57617
ORIGIN SPECIES	Human
ACCESSIONS	BC001513
VECTOR	pDONR223
NCBI SYMBOL	VPS18
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]