

GNA11-guanine nucleotide binding protein (G protein), alpha 11 (Gq class) Gene

Reference: PTXBC001528

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FBH; FBH2; FHH2; GNA-11; HHC2; HYPOC2; guanine nucleotide-binding protein subunit alpha-11; g alpha-11; guanine nucleotide binding protein (G protein), alpha 11 (Gq class); guanine nucleotide-binding protein G(y) subunit alpha; G protein subunit alpha 11		

For research use only.

GENE DESCRIPTION	guanine nucleotide binding protein (G protein), alpha 11 (Gq class)
GENE ID	2767
ORIGIN SPECIES	Human
ACCESSIONS	BC001528
VECTOR	pDONR223
NCBI SYMBOL	GNA11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatcgctgtgctgcttgcgccgcccacacggggaccctgcacggctgcttctggcctcgacagatgacaaaagaacagccccaaatacgaccactccaaccagcagttcccgcctgctgcccgcactgtcaggcctgcccggcctcctcgtccgagggtgtctgc
tggcctctggggcagaagagcggggagccccgtggaagggtcaggggagaccaggtcaggcgagctacattctggtgatcagccccatggggagacgggctggcgggataccgccccggccttccccacacacttctgtctacccgggaagcgtctttttgtgccaggt
gtctacctaagaggggtgtgccagaagcccccatggcgagtgctggggcccgcggtgcccctgggggagcagatggggccacccctggcagggcgctacaactttttcagcagcgagccctctggggggcctgtgcttgggcactctga