

GDE1-glycerophosphodiester phosphodiesterase 1 Gene

Reference: PTXBC001538

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	363E6.2; MIR16; glycerophosphodiester phosphodiesterase 1; RGS16-interacting membrane protein; membrane interacting protein of RGS16		

For research use only.

GENE DESCRIPTION	glycerophosphodiester phosphodiesterase 1
GENE ID	51573
ORIGIN SPECIES	Human
ACCESSIONS	BC001538
VECTOR	pDONR223
NCBI SYMBOL	GDE1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtggctgtggaggaccaggcgccctcctggcccttctcctcctgctgctagtgtgctgctggtgacgagcccggtcaatgcctgcctcaccggcagcctctcgttctactgcgcgtcttcagctttgagccggtgccctcttcagggccctgcaggtgtcaagccc
cgggaccgcatttgcctatgccaccgtggcgagccagcagcgcggcggagacacgctggcgccattcggcagctaagaatggagcaacaggcgtggagttggacattgagtttactctgacgggattctctcttaa