

HOXC11-homeobox C11 Gene

Reference: PTXBC001543

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HOX3H; homeobox protein Hox-C11; homeo box C11; homeobox protein Hox-3H; homeobox C11		

For research use only.

GENE DESCRIPTION	homeobox C11
GENE ID	3227
ORIGIN SPECIES	Human
ACCESSIONS	BC001543
VECTOR	pDONR223
NCBI SYMBOL	HOXC11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtttaactcggfcaacctgggcaactctctcgcgcgaaggagggggcgagatttcggcgagcgaggagctgcgcctcaacctctatctgccagttgcacttactacatgcccgagttctcaggttctctctctgccccaggccccctctcgtcagatctctac
cctactcggcccaagtgcctccgggaggtctctacggcctggagccatccggcaagtgccaccatcggaacagctactctctctgctatgcggcgccgagcagcttatgcaccgggagtgctctctctccacgtcaccgagatctctatgaaaaacgaagctctac
ggcggccaccaccacccagcggcccgacgcaacccccgggcttctactctcagfcaacaagaacagcgtctgctcaagcctcgaccgtttcttgacaacgctactgcggfggcgccgacccgccgagccccctgctccggcaaggcgagggccaaggggg
agcccaggcaccgccctcgggactggcgtccgggctgaggcggtgcgagggcgaggctgaggaggagaacaaatccagctcgtccggttcagcccactcgtggccaaggagccggccaaggagcggccccaacgcccccgaccgcaagaagcgt
gcccctattcgaattccagatccgggaactggagcgagagttttctcaacgtgtatatacaaaaagagaagcggctgcagctgtcccgatgctgaacctgacggaccgacaagtgaatttggfctcagaacagaaggatgaagaaaaagaactgagcagagaccggctgcagta
tttctcgggaatctctcgtctaa