

CHMP2B-chromatin modifying protein 2B Gene

Reference: PTXBC001553

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ALS17; CHMP2.5; DMT1; VPS2-2; VPS2B; charged multivesicular body protein 2b; VPS2 homolog B; chromatin modifying protein 2B; vacuolar protein-sorting-associated protein 2-2		

For research use only.

GENE DESCRIPTION	chromatin modifying protein 2B
GENE ID	25978
ORIGIN SPECIES	Human
ACCESSIONS	BC001553
VECTOR	pDONR223
NCBI SYMBOL	CHMP2B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgctcccttcaagaagaaacgctggatgatgtaataaaggacagatcgagagttacgaggtacacagagggctataatcagagatcgagcgtttagagaacaagaaacagctggaattagaattaagaaatggccaagattggaataaggaagcttgcgaagtttagccaacaacttgcctctacggaaacagaagacgagaacttttctgtaagtcaaaagttaacttctatgtctacacaaacaaagtgaatgccaaatgaagatggctggagcaatgtctaccacagcaaaacaatgcaggcagttacaagaagatggatccacaaagacattacaacaatgcagaattccagaaggaacatgaaatggaatgactgaagaaatgatcaatgatacacttgatgacatctttgacggttctgatgacgaagaagaaagccagatattgtgaatcaagtctctgatgaattggaattgaaattctggaagatggccaagctccatcagctgctcgaagcttaccatctgcctctacttcaaaggctacaatctcagatgaagagattgaacggcaactcaaggctttaggagtagattag