

BNIP3L-BCL2/adenovirus E1B 19kDa interacting protein 3-like Gene

Reference: PTXBC001559

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BNIP3a; NIX; BCL2/adenovirus E1B 19 kDa protein-interacting protein 3-like; BCL2/adenovirus E1B 19 kDa protein-interacting protein 3A; BCL2/adenovirus E1B 19-kd protein-interacting protein 3a; BCL2/adenovirus E1B 19kDa interacting protein 3 like; NIP3-lik		

For research use only.

GENE DESCRIPTION	BCL2/adenovirus E1B 19kDa interacting protein 3-like
GENE ID	665
ORIGIN SPECIES	Human
ACCESSIONS	BC001559
VECTOR	pDONR223
NCBI SYMBOL	BNIP3L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctgccacactagtcgagccgcccgccttgcacaacaacaacaactgcgaggaaatgagcagttcttgcctccgcccgcctcaacagttctgggtggagctacctgaacagcagcaatggcaatgataatggcaatgggaaaaatgggggctggaacac
gtaccatcctatcctccatccacaatggagacatggagaaatcttttggatgcacaacatgaatcaggacagagtagttccagaggcagttctcactgtgacagcccttcgccacaagaatgggcagatcatgtttgatgtggaatgcacaccagcaggaccatagctctcagtc
gaagaagaagtgtgtagaggagagaaggagtcgagcctttgaagaaaagtcgagctgggtatcagactgtccagtagacccgaaaacattccaccaaggagttccacttcagacaccctaacgttctgtgtcttaagcatgaggaaaagtgagccatgaagaagggtatt
ttctccgcagaattctgaaggtgttcattccatctctcttctcatgttttggcttgggctaggcatctatattgaaagcgactgagcacacctctgcagacactactga