

APTX-aprataxin Gene

Reference: PTXBC001628

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AOA; AOA1; AXA1; EAOH; EOHA; FHA-HIT; aprataxin; forkhead-associated domain histidine triad-like protein		

For research use only.

GENE DESCRIPTION	aprataxin
GENE ID	54840
ORIGIN SPECIES	Human
ACCESSIONS	BC001628
VECTOR	pDONR223
NCBI SYMBOL	APTX
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaggaccccaaatgcaggtttacaagatgagcaggtgggtgataaaggataatacccaaggcccggtaccattggctgtcttaccgtggacaccttcattccagtctgaaggctgtggccagggaacacctgaactccttaagcatatgcacactgtggggaaaaggat
tgtagattttgctgggtccagcaaaactccgctccgattgggtaccacgcccattccgagatgagccatgtacatctcatgtgatcagccaggattttgattctccttgcctaaaaacaaaaacattggaattctttcaatacagaatacttctagaatcacaagctgtgatcagatggtaca
agaggctggtagagtaactgtccgagatgggatgctgagctcttgaagctgccccttctgtcatgagtgccagcagctgctgccttcattcctcagctgaaagaacatctcaggagcactggacacagtgga