

RPA2-replication protein A2, 32kDa Gene

Reference: PTXBC001630

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	REPA2; RP-A p32; RP-A p34; RPA32; replication protein A 32 kDa subunit; RF-A protein 2; replication factor A protein 2; replication protein A 34 kDa subunit; replication protein A2, 32kDa; replication protein A2		

For research use only.

GENE DESCRIPTION	replication protein A2, 32kDa
GENE ID	6118
ORIGIN SPECIES	Human
ACCESSIONS	BC001630
VECTOR	pDONR223
NCBI SYMBOL	RPA2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtggaacagtggattcgaagctatggcagctcctacacgggggagcggcgctacacgagtcacccggggggttggatgccccgaccttcaagccgaaaagaatcaagagccgagccagcacattgtgccctgtactatatctcagctgttttgcactttgggtga
tgaagtgttcagaattgggaatgttgagatttcacaggtcactatttggggatcatcagacatgcagagaaggtcccaaccaacattgtttacaaaatagatgacatgacagctgcacccatggacgttcgccagtggttgacacagatgacaccagcagtgaaaacactgtgttctcctca
gaaacatatgtgaaatggcagggccacctgagatctttcagaacaaaaagacctggtagccttfaagatcatgccccggaggatgatgaatgagttcaccacacatatctggaagtgatcaatgcacacatggttactaagcaaaagccaacagccctcagcaggagagcaccta
tcagcaatccagggaatgagtgagcaggggaacttgggtgggaatagcttcatgcccaaatggcctcactgtggcccaaacacaggtgttgaaattgattaaggctgtccaagacctgaagggtgaacttcaggatcaagaaccagctgaaacacatgtctgtatcctcaatcaagca
agctgtgattttctgagcaatgaggggcacatctattctactgtggaatgatgacatttaaatccacagatgcagaataa