

NDUFV2-NADH dehydrogenase (ubiquinone) flavoprotein 2, 24kDa Gene

Reference: PTXBC001632

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CI-24k; NADH dehydrogenase [ubiquinone] flavoprotein 2, mitochondrial; NADH dehydrogenase (ubiquinone) flavoprotein 2, 24kDa; NADH dehydrogenase ubiquinone flavoprotein 2, mitochondrial; NADH-ubiquinone oxidoreductase 24 kDa subunit; NADH-ubiquinone oxido		

For research use only.

GENE DESCRIPTION	NADH dehydrogenase (ubiquinone) flavoprotein 2, 24kDa
GENE ID	4729
ORIGIN SPECIES	Human
ACCESSIONS	BC001632
VECTOR	pDONR223
NCBI SYMBOL	NDUFV2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttctctcgcggcgctccgggcccggcgctgctcaccgccactggggaagacatgaaggaaattgcataagacagttatgcaaaatggagctggaggagctttattgtgcacagagatactcctgagaataacctgatactccatttgattcacaccagaaaactataaga
ggatagaggcaattgtaaaaactatccagaaggccataaagcagcagctgttctccagtcctggatttagcccaaaggcagaatgggtgggtggccatctctctatgaacaaagggtgcagaagtgtttacaagtagcctccaatgagagtatatgaagtagcaactttttatacaatgtataatc
gaaagccagttggaagtatcacattcaggtctgcaactactacacctgcatgcttgaactctgacagcactatggaggccattcagaaaaagcttgaataaagggtgggagactacacctgacaaacttttactcttatagaagtggaaatgtttaggggcctgtgtgaacgaccaat
ggttcaataaatgacaattactatgaggattgacagctaaggatattgaagaattatgatgagctcaaggctggcaaatccaaaccaggcccaaggagtgagcgtcttctgtgagccagctggagggtcttacctttgactgaaccaccaaggacctggatttggtgtacaa
gcagcctttaa