

VAMP8-vesicle-associated membrane protein 8 (endobrevin) Gene

Reference: PTXBC001634

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EDB; VAMP-8; vesicle-associated membrane protein 8; vesicle-associated membrane protein 8 (endobrevin); vesicle associated membrane protein 8		

For research use only.

GENE DESCRIPTION	vesicle-associated membrane protein 8 (endobrevin)
GENE ID	8673
ORIGIN SPECIES	Human
ACCESSIONS	BC001634
VECTOR	pDONR223
NCBI SYMBOL	VAMP8
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggagccagtgagggtggaggaatgacgtgtgctgggaacctgcaagtgagggtggaggaggttaagaatattatgaccagaatgtggagcggtatctggccgggggggaaaacttgaacatctccgcaacaagacagaggatctggaagccacatctgagcacttcaag
acgacatcgagaaggtggctcggaattctgtggaagaacgtgaagatgattgtccttatctgctgattgttttatcatcatcctctcattgtgctctttgccactggtgccttctcttaa