



Reference: PTXBC001661

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ABC28; EST133090; HUSSY-18; HUSSY18; ATP-binding cassette sub-family F member 2; ABC-type transport protein; ATP-binding cassette, sub-family F (GCN20), member 2; iron-inhibited ABC transporter 2; ATP binding cassette subfamily F member 2		

GENE DESCRIPTION	ATP-binding cassette, sub-family F (GCN20), member 2
GENE ID	10061
ORIGIN SPECIES	Human
ACCESSIONS	BC001661
VECTOR	pDONR223
NCBI SYMBOL	ABCF2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]