

## FLI1-Friend leukemia virus integration 1 Gene

Reference: PTXBC001670

|                |                                                                                                                                                                                                                     |                     |                |
|----------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                 | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | EWSR2; SIC-1; Friend leukemia integration 1 transcription factor; Ewing sarcoma breakpoint region 2; Friend leukemia virus integration 1; transcription factor ERGB; Fli-1 proto-oncogene, ETS transcription factor |                     |                |

**For research use only.**

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | Friend leukemia virus integration 1                                                              |
| <b>GENE ID</b>          | 2313                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC001670                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | FLI1                                                                                             |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atggacgggactattaaggagctctgtcgggtgagcagcagaccatgcccttttgactcagcgtacggagcggcagcccatctcccaaggccgacatgactgcctcggggagtcctgactacgggcagccccaagaatcaacccctccaccacagcaggagtgatcaatc  
agccagtgggggtcaacgtcaagcgggagtgatgaccacatgaatggatccagggagtcctcgggtgacgtcagcgttagcaaatgcagcaagctgtgggcggaggcagtcacccccatgaactacaacagctatggacgagaagaatggccccctctcccaacatgacca  
ccaacgagaggagagtcacgtcctcccgagacccccacactgtggacacagagcatgtgaggaatggcgtggagtcggccataaaggagtagacgttggatggagatcgacacatccttttccagaacatggatggcaaggaaatgtgtaaaatgaacaaggaggacttctccgcgc  
caccacctctacaacacgggaagtgtgtgtcacacctcagttacctcagggaagttaactgtcgtgcctataatacaacctccacaccgaccaatctcacgattgagtgtaagaagacccctcttatgactcagtcagaagaggagctggggcaataacatgaattctggcctcaac  
aaaagtctcccttggaggggcacaaacgatcagtaagaatacagagcaacggccccagccagatccgtatcagatcctggggccgaccagtcgctagccaacctggaaagcgggcagatccagctgtggcaattcctcctggagctgctctccgacagcggcaacggcagc  
tgtatcacctgggaggggaccaacggggaggtcaaatgacggccccgatgagtggtggccagcgcgtgggcgagcggaaaaagcaagcccaacatgaattacgacaagctgagccggccctcgttattactatgataaaaacattatgacaaagtcacggcaaaagatgctt  
acaaatttgactccacggcattgccagggtctgcagccacatccgacggagtcgtccatgtacaagtaccctcttgacatctctacatgccttctaccatgccaccagcagaaggatgaaattgtccctcccatccatccatgctgctcctccagctcttctggagccgcatc  
acaatactggacctccccacgggggaattaccccaacccaacgtcccccgcacatcaacacccacgtgcttcacacttaggcagctactactag