



Reference: PTXBC001676

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IFNRP; SKMc15; SM15; interferon-related developmental regulator 2; interferon-related protein; interferon related developmental regulator 2		

For research use only.

GENE DESCRIPTION	interferon-related developmental regulator 2
GENE ID	7866
ORIGIN SPECIES	Human
ACCESSIONS	BC001676
VECTOR	pDONR223
NCBI SYMBOL	IFRD2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]