

C19orf42-chromosome 19 open reading frame 42 Gene

Reference: PTXBC001680

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0608 protein C19orf42; C19orf42; small integral membrane protein 7		

For research use only.

GENE DESCRIPTION	chromosome 19 open reading frame 42
GENE ID	79086
ORIGIN SPECIES	Human
ACCESSIONS	BC001680
VECTOR	pDONR223
NCBI SYMBOL	C19orf42
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccgagtcctcattcctcctgatgactcaagaatgttttaccagaaaaccgacaacctccagaaaagtcgaagctcgtgtgggtgggaaaagtgttcgccgaggtgtgcatggttccagccacgtccctgtttcaaatagatag