



Reference: PTXBC001689

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CAC; CACT; mitochondrial carnitine/acylcarnitine carrier protein; solute carrier family 25 (carnitine/acylcarnitine translocase), member 20; solute carrier family 25 member 20		

For research use only.

GENE DESCRIPTION	solute carrier family 25 (carnitine/acylcarnitine translocase), member 20
GENE ID	788
ORIGIN SPECIES	Human
ACCESSIONS	BC001689
VECTOR	pDONR223
NCBI SYMBOL	SLC25A20
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]