

C11orf59-chromosome 11 open reading frame 59 Gene

Reference: PTXBC001706

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RhoA activator C11orf59; C11orf59; PDRO; Ragulator1; p18; p27RF-Rho; ragulator complex protein LAMTOR1; late endosomal/lysosomal adaptor and MAPK and MTOR activator 1; lipid raft adaptor protein p18; p27Kip1-releasing factor from RhoA; p27kip1 releasing f		

For research use only.

GENE DESCRIPTION	chromosome 11 open reading frame 59
GENE ID	55004
ORIGIN SPECIES	Human
ACCESSIONS	BC001706
VECTOR	pDONR223
NCBI SYMBOL	C11orf59
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggggtgctgtacagcagcagagaacgaggactcggaccaggaccgagaggagcgggaagctgctgctggaccctagcagccccctacaaagctctcaatggagccgagcccaactaccacagcctgcttccgctgcactgatgagcagccctgctcttccatccttgcca
agacagccagcaacatcattgatgtgtctgtcagactcacaggcagtgagcagcatgagtacatggaccgtgccagcagtlacagcaccgctgtgctgtgctgagcagcagcctgacccattggaagaagctgccaccgctgccgtctctaccagcagcccaacgaagtgtg
gccagtgcagccatcccgcttctgtattgcagcaggtctccaggtatgctgtatgcctacagtgcaatttctcagatccgtgtggacgcaaaagaggagctggtgtacagtttgggatcccatga