

C14orf166-chromosome 14 open reading frame 166 Gene

Reference: PTXBC001722

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0568 protein C14orf166; CGI-99; CGI99; CLE; CLE7; LCRP369; RLLM1; hCLE1; CLE7 homolog; RLL motif containing 1; chromosome 14 open reading frame 166		

For research use only.

GENE DESCRIPTION	chromosome 14 open reading frame 166
GENE ID	51637
ORIGIN SPECIES	Human
ACCESSIONS	BC001722
VECTOR	pDONR223
NCBI SYMBOL	C14orf166
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttcgacgcaagttgacgctctgactaccacaaccccgccgcttcaactgcaaagatgaacagaatttagaaactcatgtttgcttgaagaccagaaaatcaggcactacaagattgaagacagaggggaatttaagaacatccacagcagcactggcccaagttcttga
aaagtatctcagagatgttaactgtccttcaagattcaagatcgacaagaagctattgactggctctgtgttttagctgttagacttgaatatggagataatgctgaaaaatacaaggatttagctacgtgataattcaaaaactgctgacaatgcaactaaaaatgcagaaccattgatcaatttggat
gtaataatcctgattttaaggctggtgtgatggctttggtaacctgcttcagattcagcgtcatgatgattacctggtaattgcttaaggcaattcggattttggttcaggagcgcctgacacaggatgcagttgctaaggcaatcaacaaaaaggaggcttacctgttcttagacaacata
tcttggtttgacacaggagatgcagttcttaataagctgctcaaatctgcgattgtgacacatagaggagctcagagagctacagacaaaatcaacgaagccatagtagctgttcaggcaattattgctgatccaaagacagaccacagactgggaaaagttggaagatga