

EGLN2-egl nine homolog 2 (C. elegans) Gene

Reference: PTXBC001723

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EIT6; HIF-PH1; HIFPH1; HPH-1; HPH-3; PHD1; egl nine homolog 2; HIF-prolyl hydroxylase 1; estrogen-induced tag 6; hypoxia-inducible factor prolyl hydroxylase 1; prolyl hydroxylase domain-containing protein 1; egl-9 family hypoxia inducible factor 2		

For research use only.

GENE DESCRIPTION	egl nine homolog 2 (C. elegans)
GENE ID	112398
ORIGIN SPECIES	Human
ACCESSIONS	BC001723
VECTOR	pDONR223
NCBI SYMBOL	EGLN2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtggcgtgttaccaggaacgggctcgggtacgtaaggcaggttgacaatccccacggcgatgggcgtgcacacgtgtatctattacgtgaatcagaactgggacgttaaggtgcatggcgccctgctgcagatcttcctgagggccggcccggtgtagccaacatcgagcca
ctctttgaccgggttctcatctttctggtctgaccggcggaacccccacgaggtgaagccagcctatgccaccaggtacgccatcactgtctggtattttgatgccaaggagcgggcagcagccaagacaagatcagctagcatcaggacagaaaggtgtccaagtacgtgtacacagcc
gctacgcccacctag