



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

FAM96B-family with sequence similarity 96, member B Gene

Reference: PTXBC001733

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CGI-128; CIA2B; MIP18; mitotic spindle-associated MMXD complex subunit MIP18; MSS19-interacting protein of 18 kDa; family with sequence similarity 96 member B		

For research use only.

GENE DESCRIPTION	family with sequence similarity 96, member B
GENE ID	51647
ORIGIN SPECIES	Human
ACCESSIONS	BC001733
VECTOR	pDONR223
NCBI SYMBOL	FAM96B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]