



Reference: PTXBC001744

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DNA-binding protein SATB1; special AT-rich sequence binding protein 1 (binds to nuclear matrix/scaffold-associating DNA); SATB homeobox 1		

For research use only.

GENE DESCRIPTION	SATB homeobox 1
GENE ID	6304
ORIGIN SPECIES	Human
ACCESSIONS	BC001744
VECTOR	pDONR223
NCBI SYMBOL	SATB1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]